

Adjustable Backrest

Lightweight Aluminium Adjustable Bed Backrest with Multi-Position Recline
Detachable Head Cushion | Anti-Slip Base | Up to 100 kg

The RUNMEX Adjustable Backrest is designed to provide comfortable and reliable back support for individuals who need to sit up in bed for reading, eating, resting, or recovering from surgery. Built with a lightweight yet durable aluminium frame, it features a multi-position reclining mechanism that allows users to adjust the angle for maximum comfort. The breathable fabric support and detachable head cushion help reduce pressure on the back, neck, and shoulders during extended use.



Key Features

- + Lightweight & durable aluminium frame
- + Multi-position adjustable reclining backrest
- + Comfortable breathable fabric support
- + Soft detachable head pillow
- + Foldable & portable design
- + Ergonomic back support
- + Strong locking adjustment mechanism
- + Anti-slip base for stability
- + Rust & corrosion resistant frame
- + Easy to clean & maintain
- + Ideal for home, hospital & rehabilitation use

Product Highlights

- | | |
|-------------------------------|------------------------------------|
| - Lightweight Aluminium Frame | - Multi-Angle Reclining Adjustment |
| - Soft Detachable Head Pillow | - Breathable Fabric Surface |
| - Foldable & Portable | - Stable Locking Mechanism |

Adjustable Backrest

Technical Specifications

Specification	Details
Product Name	RUNMEX Adjustable Backrest
Product Type	Adjustable Bed Backrest
Brand	RUNMEX
Frame Material	Lightweight Aluminium Alloy
Frame Finish	Silver Anodized Finish
Fabric Material	Breathable Oxford/Nylon Fabric
Headrest	Soft Removable Cushion
Reclining Positions	Multiple Adjustable Angles
Adjustment Mechanism	Locking Support System
Foldable	Yes
Weight Capacity	Up to 100 kg
Product Weight	Approx. 2-3 kg
Colour	Blue & Silver
Base	Anti-Slip Support Frame
Usage	Indoor
Assembly	Ready to Use
Package Contents	1 x Backrest, 1 x Head Cushion, 1 x User Manual
Origin	Made in India - Bhajanpura, Delhi

Ideal For

- Elderly Individuals
- Patients on Bed Rest
- Post-Surgery Recovery
- Neck & Back Pain Relief
- Orthopedic Rehabilitation
- Pregnant Women
- Home Healthcare
- Hospitals & Nursing Homes

Safety Instructions

- Place the backrest on a stable mattress before use
- Ensure the locking mechanism is securely engaged before leaning
- Do not exceed the recommended weight capacity
- Check adjustment joints regularly for secure operation
- Clean the fabric with a damp cloth and mild detergent
- Store in a dry place when not in use

[Request a wholesale quote](#)

Call 8527594730 / 7042452120 | Email runmexteam@gmail.com